

Tutorial 1: First *In-silico* microscopy image

Subhamoy Mahajan

22 Dec, 2020

1. Generate the the PSF:

The point spread function (PSF) can be generated using the following command

```
term$ python run_genpsf.py
```

It will create two PSF for wavelength 670 nm and 518 nm (“PSF_gandy_lam518_fs800.dat” and “PSF_gandy_lam670_fs800.dat”). The code is currently very slow. I will work on GPU accelerations in later versions (or hopefully someone else can help me with that).

2. Generate in-silico monochrome images.

(a) generating Image data files

The image data file containing resultant fluorescence intensity for each pixel can be calculated using the following commands,

```
term$ ../../gen_mono -p parameters.dat -f dp100.gro -o img100
```



```
term$ ../../gen_mono -p parameters.dat -f dp2000.gro -o img2000
```

It will create two pair of files “img100_lam670_fs800.dat”, “img100_lam518_fs800.dat”, “img2000_lam670_fs800.dat”, and “img2000_lam518_fs800.dat”.



(b) Render monochrome images

Monochrome images can be rendered using the following commands,

```
term$ python ../../render_mono.py -f img -p png_param.dat -t 100
```




```
term$ python ../../render_mono.py -f img -p png_param.dat -t 2000
```



3. Generate colored in-silico monochrome image.

```
term$ python ../../mono2color.png -f img -p png_param.dat -t 100
```



```
term$ python ../../mono2color.png -f img -p png_param.dat -t 2000
```

