

# Tutorial 1: First *in-silico* microscopy image

Subhamoy Mahajan

22 Dec, 2020

## 1. Generate the the PSF

The point spread function (PSF) is generted using the following command

```
term$ python run_genpsf.py
```

It will create two PSF for wavelength 670 nm (“img100\_lam670\_fs800.dat”) and 518 nm (“img100\_lam670\_fs800.dat”).

## 2. Generate *in-silico* monochrome images.

### (a) Image data files

The image data file containing resultant fluorescence intensity for each pixel can be calculated using the following commands,

```
term$ ../../gen_mono -p parameters.dat -f dp100.gro -o img100
```


```
term$ ../../gen_mono -p parameters.dat -f dp2000.gro -o img2000
```

It will generate two pairs of files “img100\_lam670\_fs800.dat”, “img100\_lam518\_fs800.dat”, “img2000\_lam670\_fs800.dat”, and “img2000\_lam518\_fs800.dat”.


### (b) Render grey-scale images

*In-silico* monochrome images can be rendeted using the following commands,

```
term$ python ../../render_mono.py -f img -p png_param.dat -t 100
```



```
term$ python ../../render_mono.py -f img -p png_param.dat -t 2000
```




It generates four files “mono\_img100\_lam670\_fs800\_I0.13.png”, “mono\_img100\_lam518\_fs800\_I0.25.png”, “mono\_img2000\_lam670\_fs800\_I0.13.png”, and “mono\_img2000\_lam518\_fs800\_I0.25.png”


### 3. Generate colored *in-silico* microscopy image.

Colored *in-silico* microscopy images can be generated using the following commands,

```
term$ python ../../mono2color.py -f img -p png_param.dat -t 100
```



```
term$ python ../../mono2color.py -f img -p png_param.dat -t 2000
```



It generates two files “img100\_fs800\_T1\_I\_0.13\_0.25.png”, and “img2000\_fs800\_T1\_I\_0.13\_0.25.png”.